

bellavista Trolley

Working with more convenience

bellavista®
the art of ventilation

The bellavista trolley (301.105.000) was designed to make working with bellavista more convenient. It is height adjustable and has a full surround mounting rail for accessories. A silicon handle allows easy storage of cables and hoses and with the larger wheels improves handling. Many optional accessories can be positioned in front or back of the trolley and are infinitely variable.

Features:

- High-quality design
- Improved stability
- Bright colours
- Very easy to assembly
- Large (125 mm) double and easy roll wheels
- Comfortable handling
- Seamless height adjustability (95 – 110 cm)
- Ergonomic, silicon coated handle
- Surround side rail to mount support arm or accessories
- Height adjustable rails to mount e.g. humidifier/filter heaters on vertical column

Additional accessories for the bellavista trolley:

Trolley Basket (302.520.000)

For your convenience, accessories can be stored at your workplace.

Oxygen Cylinder Holder (301.551.000)

Our Oxygen cylinder holder can be used for two cylinders each up to 2l. Convenient for intra-hospital transports.

Humidifier Bracket Single (304.071.000)


For mounting a humidifier to the Trolley.

Humidifier Bracket Double (304.072.000)

For mounting of two devices e.g. humidifier and filter heater.

Mounting Rail (304.073.000)

Additional rail mount for devices with a mounting clamp.



bellavista Support Arm

Working with more convenience

bellavista®
the art of ventilation

High-end support arm (301.106.000) for bellavista patient circuits.

Features:

- High-end Swiss quality design
- Single joint fixation
- Longer reach
- Side rail mounting
- Sturdy fixation

The support arm includes:

- Dual Limb circuit holder
- Rail mount clamp



Additional accessories for the support arm:

Single Tube Holder (302.790.000)

Special for single limb and coaxial circuits.

Tube Support Holder (302.792.000)

additional support for parts like flow sensor or capnography adapters.



imtmedical

imtmedical ag . Gewerbestrasse 8 . 9470 Buchs . Switzerland
T +41 81 750 66 99 . www.imtmedical.com

All ® and ™ labeled brands and products are property of the company imtmedical and can not be used without written consent. Subject to changes without notice.